

New & Refurbished

# Veterinary Equipment

Supplies, Services & Support in Canada



*Medna*   
Medical Equipment Inc.

Offices BC & ON: 833 333 2581  
[sales@mednamedical.com](mailto:sales@mednamedical.com)  
[andrew@mednamedical.com](mailto:andrew@mednamedical.com)

# Top 4 Reasons Veterinarians Across Canada Choose **Medna** Medical Equipment.

1. The “Medna Customer Bill of Rights” ([See Page 32 Now for Full Details!](#))
2. Competitive Pricing, Expert Advice, Friendly Ongoing After Sale Support.
3. Medna is the Canadian Service Depot for 100% of the Products Represented  
(Avoid the Unnecessary, Frustrating, Service, Repair/Warranty Runaround).
4. No Cross-Border Nightmares... No Unexpected Customs Delays; No Costly Hidden Brokerage Fees  
([Ask us for a cost breakdown comparison direct from USA](#)).

Medna Represents The Top  
Veterinary Equipment  
Manufacturers & Suppliers  
Across Canada

---

## More Important Reasons Leading Veterinarians Across Canada Choose Medna

- Same Day Shipping even large capital items like X-Ray Dental Units, etc.,
  - You start earning income from what you just paid for fast.
  - Local Support - Offices in BC, Ontario Alberta and Ottawa.
- 

## Did You Know Medna Specializes In Everything X-Ray?

- New Clinic Design, Specifications, Supply, Installation, Training, Support.
- HARP, RPS Surveys, Calibration, Repair.
- Digital Flat Panels, Service Repair Plus Maintenance Service Agreements.
- Canadian Compliant, Secure DICOM Cloud Backup System & Support.

For fast Local Support, Expert Equipment Advice, Local Service or Maintenance  
**Call Toll Free 833 333 2581 or Request Quote Online at [MENAMEDICAL.COM](http://MENAMEDICAL.COM)**

# MEDNA

## VETERINARY

# Contents

2	X-ray Systems
3	X-ray Solutions DR/CR
4,5	CT
6	C-Arm
7	X-ray Certificate
8	X-ray Accessories
9	Dental X-ray
10	Dental Software
11	Dental Care
12	Oxygen Concentrators
13	Lights
14	Monitors
15, 16	Services
17	Ultrasound
18	Pumps
19	Autoclaves
20	Warmers
21	Instruments
22	Supplies & Accessories
23,24	Run Doors, Cages & Kennels
25, 26, 27, 28	Tables: Exam & Surgery
29	Bill of Rights
30	Special Offers!

# Veterinary X-Ray System

Software  
Includes TTA  
VHS, TPLO  
Measurements

View your images anywhere, anytime,  
on any internet enabled device

- . Web PACS: plus 6 months unlimited cloud Storage and user access included.
- . Remotely View, store, send and receive your DICOM images 24/7/365.
- . Use a clinic iPad or tablet education & diagnosis review.

## 4-Way Float Top Table

- . Dimension:  
64 x 28"
- . Optional:  
SMART collimator with ATD (Automatic Thickness Detection)  
Easy clean-up surface with urine trap

Optional Instant  
Preview Capability



## Flat Panel Detector (5 Year Warranty)

- . 14x17 & 17x17 tethered or wireless
- . DIRECT Deposit (Csl AED) for  
low dose, high-quality diagnostic  
images, every time

## Easy Touch Intuitive Control Panel

- . Preset Factors: Or use the  
one-tap customization
- . Optional Image Preview:  
Stay table side and save time



# X-Ray Components

High - Voltage Cable



Collimator



X-Ray Tube



High Frequency  
Generator



\* These items can be sold separately or as a complete unit




# Digital X-Ray Solutions

Bi-Lateral  
Communication  
Software

## DR Flat Panel & CR Scanner

- . Superior image quality captured in 2-3 seconds
- . Viewing & Measurement Tools: Zoom, contrast, brightness, side-by-side, etc.
- . Measurements: TTA, VHS, TPLO
- . Full DICOM compliance with PACS software
- . Capable of transmitting images via e-mail and IP
- . Archival long-term storage and retrieval capabilities
- . On site & remote digital application training
- . Instant delivery
- . Exemplary technical support in Canada



**Upgrade to Digital on  
Any Existing X-Ray System  
Regardless of its Age  
In Only Half a Day!**

### OPTIONS INCLUDE

- . TETHERD
- . WIPELESS
- . 10x12
- . 14x17
- . 17x17

**Low Dose Direct  
Deposit Csl, Second  
Generation AED with  
5 Years Warranty**



mini-PC



Equine Package



CR Scanner

- . Standalone, Industrial-grade computer
- . HIPPA - Compliant encryption
- . ARM - based security module
- . Medical grade work stations and monitors available

**Go Digital For Only \$ 1.00  
No Payments up to 6 Months!**

## Software for Image Acquisition, Processing, Diagnostics & PACS Storage

- . CaptureX imaging software with diagnostics.
- . OR technology for optimized workflow.
- . 20/20 Imaging Opal RAD
- . Nexus Varex - Nexus DR
- . DIRA Veterinary

# NewTom 5G XL

## ADVANCED 3D VETERINARY IMAGING

### Outstanding quality and innovation,

### exceptional features.

NewTom 5G XL, a CBCT device incorporating unique technology, guarantees advanced top-level diagnostics by maximising patient stability and providing a broad selection of FOVs, resulting in 3D and 2D images of the very highest quality.

CBCT technology combines low X-ray doses with high spatial resolution: the outstanding quality of 5G XL images allows veterinary specialists to identify pathologies and fractures that would otherwise be unidentifiable. Patient health - enhanced by the ECO Dose mode and exclusive SafeBeam™ technology - is a key design priority.

The native 21 x 19 cm FOV and numerous smaller high resolution FOVs constitute a toolbox that is perfect for veterinary diagnostics, which demand sharp and extremely detailed images.



Ask Us  
about our  
Refurbished  
Inventory



**SUPERIOR  
DIAGNOSTIC  
QUALITY**



**MOTOR-DRIVEN  
PATIENT TABLE**



**MINIMUM TIMES,  
MAXIMUM  
USER-FRIENDLINESS**



**VETERINARY  
SOFTWARE**


## SUPERIOR DIAGNOSTIC QUALITY

Superior standard  
3D examinations  
with a device  
designed for excellent  
performance.

Ultra-high definition volumetric images with native isotropic voxel resolution and minimal artifact presence. Thanks to CBCT technology, 5G XL offers fast scanning and extremely contained X-ray doses, increasing safety for both operator and patient and ensuring a more efficient workflow.

The outstanding versatility of 5G XL makes it the perfect device for otorhinolaryngological and muscular-skeletal diagnoses, with or without radiocontrast agents, for complete upper airway analysis and for the investigation of dental and maxillofacial pathologies.

### COMBINING TECHNOLOGY, PERFORMANCE AND SAFETY



5G XL offers outstanding results with low radiation doses for the patient, all thanks to undeniable constructive excellence:

The X-ray generator, the only one of its kind in the world, with rotating anode and a 0.3 mm focal spot, optimises performance by adapting emissions to the specific needs of the examination. It also provides high filtration levels to protect the animal from more harmful radiation types.

The large Flat Panel sensor with high signal-noise ratio improves final image quality, while collimation limits exposure to the specific regions of interest.


Innovative volumetric reconstruction algorithms allow the user to control the image chain and maximise diagnostic potential while minimising the presence of artifacts.

## MOTOR-DRIVEN PATIENT TABLE.

The NewTom 5G XL is characterised by a motor-driven, carbon fibre patient table. Controlled via a panel on the machine, it ensures comfortable, easy access to the gantry for operator and patient alike.



Motor-driven table control is even easier when performed via the PC as it allows the user to adjust patient positioning directly from the workstation. The open gantry makes access to the scan area easy and simplifies, along with the three reference lasers, the management of initial patient positioning (dorsal, ventral, rostral, cranio-caudal or caudo-cranial).



360° reconstruction



eXtra FOV vision




HiRes analysis

# Skan C VET

## High Frequency Veterinary Mobile Surgical C-Arm



Ask  
Us



Mobile, High-Frequency Surgical C-Arm with 2.2kW Monoblock, 9" Tri-Mode Image Intensifier, kV Range of 40-110kV, 10.4" Touch-Screen Console with Image Adjustment Features, On-Board Storage For Up to 100,000 Images, 1k x 1k High Resolution Camera, Road Mapping & DSA (Standard), Laser Aimer, DICOM 3.0 with PACs integration, Dual 18" High Contrast Monitors and Sleek On Board Keyboard.

### Radiation:

- \*Safest C-arm in its class with near zero leakage radiation
  - \*Lower dosage with excellent image quality
- Generator:
- \*High-performance generator with up to 250kHz switching frequency
  - \*Customized and in-house developed X-ray tube

### Modes:


- \*Pulse mode for a lower dosage
  - \*DSA & Roadmap comes as a standard feature
- Service support:
- \*Best-in-the-industry warranty terms
  - \*Remote diagnostics capability
- Low dosage patient positioning:
- \*Virtual collimation for ROI selection w/o X-ray
  - \*Laser Aimer for initial patient positioning
- Workstation capability:
- \*More than 100k images storage
  - \*DICOM 3 compatible; PACS integration
  - \*2-hour inbuilt power back-up
  - \*Upgradable option for wifi capability for PACS
  - \*Upgradable option for external monitor



# + Web PACS CaptureX Unlimited Viewers

## Reasons Why Veterinary Clinics Love Using The Medna Web PACS Solution:

- . Easy to Use
- . Automatic Backups!
- . Saves You & Your Staff Time and Money
- . Access your Studies in Seconds - From Anywhere!
- . Unlimited Viewing Stations
- . Works on Any Device (PC, Mac, iPhone, iPad, Android, Tablets etc...)



## X-Ray Certificate

### Mandatory Collage Requirement

- . Radiation Protection Survey (RPS) & HARP
- . Radiation Safety Certificate
- . Certificate of Compliance



**Go to Mednamedical.com**  
and book your appointment to  
receive 10% off

\* Please note: Other names are used for this test. Please contact us for details.



# Veterinary Portable & Equine X-Ray

## PXP20HF Plus & VET-20BT

- . 40-100 KV @ 1.0 kv increment
- . Lightweight
- . Small size
- . Freedom of use
- . Zero clutter
- . Point and shoot



Lightest,  
smartest and  
most durable  
portable  
X-Ray



# X-Ray Accessories

Leaded window frames



Lead Sheetting



Lead glasses



Lead Apron & Glove Rack



Lead Thyroid Collar



Lead Apron with Integrated  
Thyroid Collar



Veterinary Lead Apparel  
and Accessories Pack



Lead Mittens



Lead Gloves (One Pair)



# Dental X-Ray

Affordable & High-Quality DC Dental X-Ray Systems



## DCX7000 Series Dental X-ray (Wall-mount & Mobile)

Introducing a power pack and intelligent HF dental system

Advanced dual microprocessor system + High Frequency DC system

- 0.4 mm focal point, Low dosage.
- Sharp and crisp image quality X-Ray.
- Fully digital-programmable operation.
- Perfect clinical and patient experience.

### Flexible Exposure System

- Triggers from distance with no-shift image position.
- Ease of movement and flexible.
- Availability of three scissors arms sizes.

### 30 Pre-programmed techniques

- Enhancing doctor's experience and flexibility.
- Exclusive LCD display.
- Streamline setting for individual.

### WE BELIEVE IN GREEN so as our system

- Clean power withdrawal while on stand by.
- Energy consumption negligible.
- Energy conservation with shorter exposure times

## DCX9000 Series Dental X-ray (Portable)

Redefining the whole portable experience with ease of movement, Safe, Lightweight, Quick transport, Ergonomic design & Convenient to use

### Single unit with Multi features

Longest Battery life in the market, Up to 800 Exposures

- 0.4 mm focal point, Super sharp output
- Dual radiation protection with Max weight of 3.0 Kg.
- Internal shielding & External backscatter shield.
- Fully digital-programmable operation.
- Fast exposure system, X-Ray field of 60mm diameter.
- FDA, CE and Health Canada Approved.

### REDEFINING INTRAORAL SENSOR

FHDX1000 (size #1), FHDX2000 (size #2) & FHDX4000 (size #4)

- Specifically designed for veterinary use
- Fiber optic plate with over 2m cable length.
- Waterproof, can be immersed in disinfectant liquid.
- Sleek CMOS sensor adjustable for every patient comfort.
- High Resolution + High Sensitivity = Accurate Diagnosis.
- Cutting Edge Image processing.
- Optimum and crystal clear image quality.
- Powerful imaging software, sec image transfer.
- Enhanced image quality, fully DICOM, Export with the option of cloud-view and storage.
- lifetime software support (manufacturer limited)



Mobile Unit



Wall Mount



Portable X-ray



FHDX2000 (size #2)




FHDX1000 (size #1)

Buy a Wall Mount or Mobile Unit with the Digital Intraoral Sensor and

**SAVE \$500**

Why  
Choose  
DC?



DS units produce a more consistent x-ray output which allows for shorter exposure times and take quality images consistently, compared to AC machines.

# Dental Software captureX Dental Software

## Software Key Features

- . Developed in Canada
- . Cloud software
- . Intuitive and smart UI (User Interface)
- . Unlimited Storage
- . Enhanced Package
- . captureX unlimited storage
- . captureX unlimited user access
- . captureX unlimited user log-in
- . captureX unlimited viewing stations
- . captureX unlimited enhanced medical report
- . captureX surgical VHS, TTA and TPLO tools
- . captureX fully DICOM format imaging
- . captureX unlimited DICOM send SCU
- . captureX unlimited email module
- . captureX unlimited download module in DICOM/JPEG
- . captureX patient USB archiving
- . captureX advanced veterinary tools
- . First year included



## Vault & View

- . Developed in Canada
- . 365/24/7 Cloud backup
- . Patient filters & Search
- . Wide variety of editing tools
- . Online access from anywhere



## Vault & View

To Register Visit:

[www.Mednamedical.com](http://www.Mednamedical.com)

# Dental Care

## Scale-Aire

Always  
in  
Stock

Canadian  
Serviced

Support  
and  
Training

- High speed, multi-function veterinary dental unit.
- All-in-one unit lets you easily switch from polishing to scaling to cutting.
- Surgical stainless steel, heat treated tips glide smoothly across tooth surfaces.
- Optional, fiber optics system with high-speed hand piece.
- Tips gently cool and lubricate teeth, and wash away dislodged debris
- Built-in 25K ultrasonic scaler maintains optimum frequency.



## Accessories & Maintenance Kits

Looking for a dental unit Maintenance Kit? Look no further we have maintenance kits for your Sonus V, Son-Mate II, Scale-Aire & Scale-Aire Mini's.

- Choose from: Basic - Deluxe packages.
- Spray Lubricants, One-Step conditioners, Drop Oil, Water Filtration Bottles
  - Cleaning Brushes, Prophy Angels & Bur Kits.
- Not seeing what you need? Contact us we will find it for you



### Dental Scaler

- Water line with filter and quick disconnect.
- Perform curettage & deep curettage.
- Scale heavy calculus and remove stains.
- Transducer (stack) inserted in the hand piece.



### Scale -Aire Mini

- A powerful and potent tool against periodontal disease.
- Connects to a nitrogen tank easily.
- Water flows internally through the tip directly to the operative site, gently cooling and lubricating the tooth, washing away the dislodged debris.



### Son-Mate II

- From 200 RPM high torque micro-motor with forward/reverse switch.
- One foot control for both scaling/polishing modes.
- Includes a scaler hand piece with an inserted stack (transducer),
  - assorted tips and a tip wrench.



# Oxygen Concentrators



## Drive Medical / 5

- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter Oxygen Concentrator was designed with enhanced quality, durability and simplicity in mind.
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device

Delivery Rate	2.0 to 10.0 LPM
Oxygen Percentage	2-10 LPM, 87-96%
Power Consumption	120 VAC, 60Hz
Power Input	639 (avg)
Warranty Limited	3 years



## Drive Medical / 10

- The 1025 Oxygen concentrator is designed to provide optimal oxygen delivery for patients with varying flow requirements.
- Assembled in the USA and built in the same shell as the 5 Liter, the 10 Liter Concentrator is also compatible with cylinder transfill systems
- Sold separately, the Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs., to the 10L Concentrator

Delivery Rate	2.0 to 10.0 LPM
Oxygen Percentage	2-10 LPM, 87-96%
Power Consumption	120 VAC, 60Hz
Power Input	639 (avg)
Warranty Limited	3 years

Little  
Maintenance

Quiet  
&  
Reliable

Call  
Today!



# Veterinary Lighting

## MDN LED lights

- 45,000lux (4181 fc) at 1 meter.
  - Uses 55% less energy compared to the traditional light sources.
  - 3 LED modules with an 20,000-hour life.
  - Adjustable focus using center handles.
- Number of LEDs 48 to 80



LED Rated Life 50,000 hours

## MDN Outpatient LED

- 90,000 lux (7990fc) at 24" (61cm).
- The unique shape of the light head reduces shadows on the examination area
- 20,000 hours average LED life

#1 Dental &  
Treatment  
Choice



All Lights Available In Ceiling,  
Wall, Floor & Double Head



Dual Head



Ceiling Mount



Mobile



Wall Mount

# Monitors



## EPM 12V VET



## EPM 10V VET



### Reliable Quality

The ePM 12M Vet comes equipped with a high-quality shell materials, fan-free design and veterinary specific accessories. They are safe, durable, and built to last:

### High-quality shell materials

Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners.

### Fan-free design

- Keeps fur out of equipment
- Reduces cross-contamination, protecting the health of patients and staff
- No internal dust, extending the life of the monitor

### Veterinary specific accessories

- Low-trauma ECG clips maintain consistent contact
- Appropriately-sized NIBP cuffs provide accurate readings
- Durability-tested for long-lasting reliability

### Reliable Quality

The ePM 10M Vet comes equipped with a high-quality shell materials, fan-free design and veterinary specific accessories. They are safe, durable, and built to last:

### High-quality shell materials

Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners.

### Fan-free design

- Keeps fur out of equipment
- Reduces cross-contamination, protecting the health of patients and staff
- No internal dust, extending the life of the monitor

### Veterinary specific accessories

- Low-trauma ECG clips maintain consistent contact
- Appropriately-sized NIBP cuffs provide accurate readings
- Durability-tested for long-lasting reliability

## Accessories

### ECG CLIP



ECG: Non-traumatic  
crocodile clips

### NIBP CUFFS



NIBP: Variety of cuffs fits  
multiple patient sizes

*Medna*  
Medical Equipment Inc.  
**Services**

For More Details & Request a Quote:  
Go to [MEDNAMEDICAL.COM](http://MEDNAMEDICAL.COM)



## Equipment Service

- Warranty (Medna is The Canadian Service Depot  
for 100% of products sold)
- Service Agreements
- Maintenance (with or without service agreement)
- Laser Equipment Service & Repair in Canada
- Loaner Equipment
- Exchange/Upgrade Programs



## X-Ray Service

- Radiation Safety Survey & Certificate RPS
- Calibration
- Repairs
- Digital X-Ray Support
- Cloud PACS Setup
- Lead Shielding & Power Requirements

**For Everything Veterinary In Canada!**

**Call 833 333 2581**

16

# MEDNA

is The Canadian Service  
Depot for 100% of products sold



## 7 Day-a-week Support in Canada

- Emergency Service 24/7
  - (Next-day support guarantee)
- On-site & Remote Technical Support
- Repairs (Why replace it when you can fix it?)
- Retrofits/Upgrades
- (Take advantage of new technology)



## Financing

- Banking Relationships
- Leasing
- 3 or 6 months no payment plan



## Training

- Staff Training
- Technique Optimization
- Digital X-Ray Training
- Dental X-Ray Training

## Consulting

- New Clinic Design: Turn-key
- Room & Radiology Preparation
- Ergonomically friendly designs

## Installation

- New & Refurbished Equipment
- Retrofits & Upgrades

**For Everything Veterinary In Canada!**

**Call 833 333 2581**

# Ultrasound



## Vetus E7



Since its founding, Mindray Animal Medical has been dedicated to developing innovative and accessible ultrasound solutions for the animal care environment. It is this commitment to excellence that has brought you ZONE Sonography® Technology+ (ZST+), evolving industrial designs and comprehensive disinfection solutions. The Vetus E7 Ultrasound System has been thoughtfully designed to overcome the obstacles clinicians face in today's challenging healthcare environment. With a sealed user interface design and revolutionary, software-based beamformer technology, ZST+, the Vetus E7 System combines best-in-class image quality with an intuitive user experience to help ensure reliable and efficient diagnosis during the most challenging exams

## DP-50 Vet



With its streamlined features and compact design, the DP-50 Vet Ultrasound System is the ideal solution for a mobile veterinary environment. The System's powerful engine and high-speed processing provide outstanding image quality and a complete information management solution. For veterinarians seeking improved ergonomics and innovative ultrasound technologies, the DP-50 Vet is an excellent choice for providing high-quality care

## Z60 Vet



The Z60 Vet Ultrasound System is a comprehensive choice for veterinarians seeking innovative and high-performing ultrasound technology. A compact System, the 260 Vet combines unsurpassed image quality with ease-of-use, and practical functionality, while meeting the financial goals of any practice.

## M6 Vet



### Intelligent Workflow

- iTouch™: Intelligent image optimization for B, color and PW mode
- Smart Tracking: Continuously track the color flow and optimize the best color box in real time scanning
- Smart Doppler: Automatic optimize the size and position of color box and doppler gate placement to ensure optimal signal color and spectral doppler sampling

### Comprehensive Functionalities

- iScape™: real-time panoramic imaging with velocity indicator, a very useful tool to evaluate the racing horse tendon
- Free Xros™ Imaging: Anatomical M mode, it can support EF measure offline and significantly shorten the exam time
- Advanced cardiac function assessment packages including Free Xros™ CM Imaging: Anatomical Curved M mode and

### Professional Design for Veterinary Use

- Veterinary information management tool for more than 6 types of animals
- Up to 22 veterinary body marks
- Abdomen, cardiac, OB, vascular, small parts and MSK exam modes and analysis packages for different animals
- Veterinary OB table
- User-defined formulas for more applications

## M9 Vet



### Dedicated Veterinary Software

- Data information
- Exam modes
- Bodymarks
- Measurements and reports

### Advanced Features

- High quality CW imaging
- TDI and quantitative analysis for confident cardiology analysis
- Tissue Tracking quantitative analysis for advanced myocardial study
- Auto measurement including Auto PW trace, Auto EF etc

## DP-30 Vet



### Dedicated Veterinary Software

- Data information
- Exam modes
- Bodymarks
- Measurements and reports

### Advanced Features

- High quality CW imaging
- TDI and quantitative analysis for confident cardiology analysis
- Tissue Tracking quantitative analysis for advanced myocardial study
- Auto measurement including Auto PW trace, Auto EF etc.



# Pumps

## EN-V5 Infusion



- 4.3 inch color touch screen, parameter setting and editing directly
- Uses any standard IV sets and stores no less than 2000 drugs in library
- Multi-function operation with 5 infusion/inclusion modes
- Infusion rates range from 0.1 to 2000ml/h
- Last therapy for rapid infusion start
- 10 hours battery back-up with Lithium polymer battery 11.1V/2600 mAh
- IP44 waterproof

## EN-V7P Infusion



- 4.3 inch color touch screen, parameter setting and editing directly
- Uses any standard IV sets and stores no less than 30 drugs in library
- Multi-function operation with 3 modes
- Electric door and anti-free flow clip
- patented design
- Double CPU for increased safety.
- Ultrasonic air-in-line detector
- Historic record more than 5000 logs
- 9 hours battery back-up with Lithium polymer battery 11.1V/2600 mAh
- IP24 waterproof, and stackable design

## EN-S7P Syringe



- 4.3 inch color touch screen, parameter setting and editing directly
- Automatic recognition of syringe sizes (5,10,20,30 and 50/60)
- Stores no less than 30 drugs in library
- Multi-function operation with 3 modes
- KVO range from 0.01 to 5ml/h, editable
- History record more than 5000 logs
- 12 hours battery back-up with Lithium polymer battery 11.1V/2600 mAh
- IP24 waterproof, and stackable design

## VV-100 Vein Viewer

- Easy to operate
- Accurate imaging
- Three operating modes



**Basic Mode**  
Basic monochromatic mode,  
clear image



**Enhance Mode**  
Two-color image,  
enhanced appearance



**Inverse Mode**  
Inverse color image,  
prominent vein



## Tubing set

- Optimally deliver IV fluids, TPN solutions, blood and blood products with a +/- 3% accuracy.
- All control-V IV pumps are now calibrated to this tubing in channel 1.
- 103",20 drops/ml, dual injection ports, rotating luer lock.
- 50 V tubing sets/cases

## IV Stands



# Autoclaves

## Clave23+ with Compressor Class S Autoclave, Non-Vacuum



Full automatic autoclave with LCD displaying,  
independent steam generator,  
Independent temperature and pressure monitoring,  
higher precision control  
System, automatic alarm for water level, door locking  
and error code  
Indicated. \*Includes Built in Printer and Datalogger  
# of Trays: 4 trays (standard equipped), max 5 trays

<b>Chamber Size</b>	F9.72"X17.7" (23L)
<b>Overall Size</b>	18.9"(W) x 18.3"(H) x 27.6"(D) (23L)
<b>Nominal power</b>	1600W
<b>Voltage</b>	110VAC - 60Hz

## Tuttnauer Manual Autoclaves



- Automatic shut off at the end of both sterilization and dry cycles.
- A long life electro polished Stainless Steel chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.

	10M	9M
<b>Chamber Size</b>	10"x18"	9"x18"
<b>Volume</b>	6gal/23L	5gal/19L
<b>Cassette Capacity</b>	3full & 3half	2full & 2half
<b>No. of Trays</b>	4	3



## Ultrasonic Cleaning Systems

- High-efficiency 60kHz cavitation for a thorough and quieter cleaning.
- Durable vinyl clad steel housing. Easy to clean and scratch resistant.
- The 1-gallon system (CSU1) holds a standard selection of instruments or 2 half sizes cassettes. The gallon holds 3 full size cassettes.

## Sterilization Ideals



Instrument Pouch Racks



Chamber Brite Autoclave Cleaner



Ultrasonic Cleaning Tablets



Water Distiller

# Warmers & Suction Pump



## Bediet Warmer

A set of 3 different reusable, machine washable warming blankets are included making it by far the best value on the market.

Remote control

- Adjust from anywhere

.Weighs 13.8lbs

.Super Quiet Operation



## Pump structure

- Diaphragm type
- Limit negative pressure value: 0.09MPa
- Negative adjustment range: 0.02 ~ 0.09MPa
- Suction rate:  $\geq 40\text{L/min}$  (terminal:  $\geq 40\text{L/min}$ )
- Suction bottle: 2500ml x 2
- Input power: 400VA
- Fuse: FLSA 250V  $\Phi 5 \times 20$

# Instruments

## High-Quality, German Surgical Instruments

### Dental Instruments

- .Surgical Blades
- .Periotomes
- .Elevators
- .Root Tip Pick
- .Extraction Forceps
- .Needle Holders

### Surgical Packs

- .Feline & Canine Spay packs
- .Standard Veterinary Eye Pack
- .Standard Veterinary Orthopedic Pack
- .Small Animal Dental Extraction Kits
- .Rabbit Dental Extraction Kits


 Medna  
Medical Equipment Inc.

## To Serve You Even Better... We now sell consumables

Because using high-quality, long-lasting and intelligently designed instruments make completing your procedures more efficient, safe and less stressful. We can provide your clinic or hospital with precision surgical instruments for every procedure at extremely competitive pricing.

Call with any important questions or to arrange a visit so you can see the instruments first hand. We distribute over 10,000 products consisting of, High Quality Surgical Implements & consumables.

#### Balfour Retractor


• 90mm • 180mm • 300mm



#### Baby Gelpi Retractor



#### Feline | Canine Spay Pack





# Supplies & Accessories

**Refurbished items available:**  
Please contact us  
for more details!

■ ■ Welch Allyn  
Diagnostic Sets



■ ■ Thermometer



■ ■ Stethoscopes



■ ■ Wall-Mounted  
Diagnostic Sets



■ ■ Thermometer



■ ■ Tono-Vet



■ ■ Digital Touch  
screen Electro-Son



■ ■ Apnea Alert



■ ■ Frameless  
Stretcher



■ ■ Cat Grabber



■ ■ Stretcher



■ ■ Positioners



■ ■ Digital  
Platform Scale



■ ■ Cat Scale





# Run Doors & Cage Accessories

**Grill Run Door**



**Glass Run Door**



**Full Frame Run Doors with Grill**



**Full Frame Run Doors with Glass**



**Oxygen Therapy Door**



**Intensive Care Unit**



**Transport Cage with Therapy Door**



**Noise Reduction Strapping**



**Clipboards**



**Cage With Pole**



**Clipboards with Bin**



**Humidifier – Single Patient Use**



**Plastisol Removable Cage Floors**



# Cages & Kennels

## All Styles of Regal Cages and Designed with these Quality Features:



**Doors** Stainless steel door spokes are deburred smooth to eliminate sharp ends. Stainless steel, electropolished removable ambidextrous door may be remounted with hinges on the right side.  
Door extends below cage floor to help keep paper liners in place.  
Door fits snug to the cage body to eliminate areas where an animal's paws could be pinched. 1875-n. (45cm) diameter stainless steel vertical spokes are spaced 1-in. apart and welded at every intersection to crossmembers and frame.  
3125-in. (79cm) diameter stainless steel door frame and crossmembers are welded together for strength.



**Latch** Stainless Steel Latch easily disassembles for sanitizing and adjusting. 3125-in. (79cm) diameter stainless steel latch keeper latches door securely. Heavy 12-gauge stainless steel latch brackets. Latch connector bar is notched for easy opening. Heavy 12 gauge stainless steel self-latching door latch connector is designed to accommodate a padlock for greater security.



**Hinges** 3125-in. (79cm) diameter stainless steel hinge pins allow easy removal of door.  
Heavy 12-gauge stainless steel die-formed barrel hinges operate smoothly for the life of the cage



**Body** 22-gauge Type 304 stainless steel body sheets are welded into one piece to create the strongest cage available. 750 in. (1.90cm) covered rear corners on top and bottom eliminate accumulation of dirt and bacteria. 750 in. (1.90cm) radius bends on all sides and back provides a smooth, easy-to clean interior surface. Distinctive floor design with fluid retaining front lip is standard on all Regal Cages. Bottom lip overlaps lower cage when used in a Regal Cage arrangement.



**Hardware** Stainless Steel installation hardware is the strongest available. Vertical molding strips and 14-gauge bolt-on corner plates mount easily with 250-20 stainless steel truss-head machine screws.



**Frame** 18-gauge stainless steel die-formed floor support brace keeps floor rigid for animal comfort. 18-gauge stainless steel die-formed rear support keeps cage level. Hidden 16-gauge stainless steel die-formed front frame makes the front of the cage extra strong.



## Sample Configurations



Single Door Stainless Steel Cages



Double Door Stainless Steel Cages



Aristo-Cat Kitty Suites, Cat Condo



Stationary



Heated



Mobile

# Exam Tables



Classic Pedestal Exam Table with Tilt Top



Classic Pedestal Exam Table with Regal 300 Electronic Scale



Classic Pedestal Exam

## Lift Table



More models & configurations available  
Contact us for details

## Premier ExamTable



## Classic Vet Scissors Lift



## Classic Tubular Base Exam Table



## Folding Wall-Mounted Exam Table





# Prep-Procedure, Wet & Dental Table

Regal Stainless Steel  
Imperial Laminated



## Options Available

**Multi Depth Sink**  
Coated Removable  
Bathing Rod




**Recessed Sink**



**Full Sink**  
Stainless Steel Removable  
Bathing Rack



## L-Arrangement, 2station





# Surgery Operating Tables



Classic  
Electric lifting with scale



Classic  
Electric lifting



Classic  
Manual lifting

## Our Selection Include:

Flat Top or V-Style Top Heated or  
Non-Heated 46" or 58" Wide Top

## 4 Easy Steps To Choosing Your Table:



STEP  
1

V-Top or Flat-Top?

STEP  
2

46" or 58"?

STEP  
3

Heated or Non-Heated?

STEP  
4

Manual or Electric?

# Tables & Surgery Accessories

Over Instrument  
Table



Utility Tables



Over Operating  
Table



IV-Stand



Chrome Base Stool



Mayo Stands



5-Caster Instrument  
Stand



Solution Stand



Each Pull-out Tray  
is sold separately



Vet Stretcher



Sponge Receptacle



Kick Bucket



Bucket



Stainless Steel  
IV Poles



Stainless Steel  
Instrument  
Tables



Operating Table  
Casters



Head Extension for  
V-Top Table



Dual-Tilt  
Retrofit Kit  
for Classic  
& Elite  
Operating Tables



Thoracic Positioner



# Customer's Bill of Rights

**Medna Medical is proud abide by the Medna Customer's Bill of Rights for the Veterinarian Industry which Protects our customers from any negative experience they might otherwise have from other veterinarian supply companies:**

- 1 We will treat our customers the way we would like ourselves to be treated, like family!
- 2 We guarantee our customers will receive honest and ethical advice, service and pricing.
- 3 The customer's best interest always take precedence in any intranction with Medna , and the customer is always right.
- 4 We will clearly explain the pros and cons of all options presented so the customer can make a wise, informed decision.
- 5 We will deliver on 100% of our promises. If we can not follow through on a commitment, we will not promise what we can not deliver!
- 6 The customer will always receive prompt communication from their representative.
- 7 Clear expectations will always be set and agreed upon by the customer and representative, to avoid misunderstandings and unnecessary problems.
- 8 All quotes will clearly state all prices and fees up front-protecting you from surprise costs like shipping, installation, duties, brokerage fees, and price changes.
- 9 Service vehicles will be fully stocked with the necessary diagnostic and repair equipment ensuring the quickest repair and shortest service times possible.
- 10 We will protect our customers with the best warranties in our industry.
- 11 We will invest in the most advanced equipment in our industry to ensure our customers get best and fastest service possible.
- 12 Our technicians will have no less than 5 years-experience before working in the field for our customers.
- 13 We will always remember to quote on all items required to have you purchase work as you expect it to do.

30

# Special Offers

Go Digital For  
**Only \$1.00** No Payments  
for **6 Months!**

## Convert Your Existing X-Ray to Digital



**SPECIAL  
DISCOUNT**

### Are You...

Waiting more than 2-3 seconds for diagnostic quality images?  
Tolerating (or embarrassed by) poor image quality?  
Retaking X-Rays?  
Exposing yourself to more radiation than necessary....?

### Then NOW is the best time to upgrade and save thousands!

Choose From Three Special Offers Plus Discounted Lease Options:  
60 Month term, 10% Buyout  
66 Month term, \$1.00 Buyout  
Skip 3/6 months and

## START MAKING MONEY NOW!

Ask us about our Loyalty Program!

Get 1 Year Of Unlimited Web PACS and Cloud  
Storage, Unlimited Viewing Stations AND A 5-Year

“**Bumper to Bumper Warranty**”



# How Does Your Veterinary Supplier Rate?

For your convenience, we've provided this important check-list with all the questions you'll want to ask before buying the products you need.

## Company Status

Has Offices, Inventory & Staff (Not 1-Man Show or Drop Shipper)  
Stock 80% of Items in Canada (Has Same/Next Day Shipping)  
Single Vendor Capable (Full Product Line & Support, including X-Ray)

## Canadian Vs. USA Supply & Service

• Local Service & Support Technicians (No "Technician Fly in Costs and Delays)  
Canadian Warranty (No Ugly Cross-Border Nightmares!)  
Canadian Service Depot (for All Products Sold)

## Shipping & Customs

No Unknown Fees for Shipping or Customs Brokerage  
Pick any shipping company you want (Some Force Expensive ones!)

## Service & Support

Next Day Support Guarantee  
Maintenance & Service Agreements on 100% of Products Sold  
Loaner Equipment Available (Avoid Costly Downtime)

	Medna vet	Big International	Online	US Dealer
Company Status	✓ ✓ ✓	✓ ✗ ✓	✗ ✗ ✗	✓ ✗ ?
Canadian Vs. USA Supply & Service	✓	?	✗	✗
Canadian Warranty	✓ ✓ ✓	? ? ✗	✗ ✗ ✗	✗ ✗ ✗
Shipping & Customs	✓ ✓	✓ ✓	✗ ✗	✗ ✓
Service & Support	✓ ✓ ✓	? ✓ ?	✗ ✗ ✗	✗ ✗ ✗