


T-Quick Autoclave



The Only Manual Autoclave with Closed Door Drying

Post Covid-19: Chairside sterilization is no longer recommended due to aerosols!

The T-Quick autoclave provides:

- * Faster cycles that deliver dry instruments and no wicking
- * Closed door drying
- * More partial loads when needed vs. full loads



Medna
Medical Equipment Inc.

Features:






- T-Quick allows you to easily track each phase of the sterilization process
- Automatic safety shutoff prevents overheating
- Water detector with automatic shutoff
- HEPA air filter provides sterile, bacteria-free air for drying.
- Aligns with AAMI guidelines regarding sterilization cycles.

T: 833 333 2581 P: 604-283-4717 C: 647 948 4797 C: 778 881 4902

E: Info@mednamedical.com

Tuttnauer
Innovation · Legacy · Partnership

Sterilization Programs

Icon	Cycle Program	Load Type	Cycle Time (min.)*
	Unwrapped 134	Unwrapped instruments	11 Minutes
	Wrapped 134	Wrapped instruments	15 Minutes
	Unwrapped 121	Unwrapped delicate instruments	25 Minutes
	Wrapped 121	Wrapped delicate instruments	25 Minutes
	Unwrapped	Single Instruments	9 Minutes

* Cycle times are based on: 1 kg standard loads (not including drying), Warm chamber (and subject to load size)

Panel Features



- ① Cycle Status Indicator
- ② Pressure Indicator
- ③ Power Switch
- ④ Sterilization Timer
- ⑤ Water Pump Switch
- ⑥ Drying and Sterilization Timer
- ⑦ Night Mode Switch
- ⑧ Temperature Selector

Technical Specifications:

- Chamber Dimensions W x L (mm) 250x467
- Chamber Volume 23L / 6.1gal
- Standard Cassettes Capacity 3 Full
- No. of Trays (4)
- Tray Dimensions WxHxD (mm) 168x20x414
- Voltage(V) / Freq.(Hz) 230V / (50/60 Hz)
- Current(A) / Power(W) 9.6A / 2200W
- Overall Dimensions WxHxD (mm) 508x362x550



Tuttnauer
Innovation · Legacy · Partnership

Medna
Medical Equipment Inc.

T: 833 333 2581 P: 604-283-4717 C: 647 948 4797 C: 778 881 4902

E: Info@mednamedical.com